Quick Set-up Guide InPro 5000 (i) Sensors

cs	Průvodce rychlým nastavením	2
da	Lynvejledning	6
de	Kurzanleitung	10
en	Quick Setup Guide	14
es	Guía de configuración rápida	18
fi	Pika-asetusopas	22
fr	Guide de paramétrage rapide	26
hu	Gyorsbeállítási útmutató	30
it	Guida alla configurazione rapida	34
ja	クイック セットアップ ガイド	38
ko	사용 매뉴얼	42
nl	Beknopte handleiding	46
pl	Podręcznik szybkiej konfiguracji	50
pt	Guia de Configuração Rápida	54
ru	Руководство по быстрой настройке	58
sv	Snabbguide	62
th	คู่มือการตั้งค่าอย่างรวดเร็ว	66
zh	快速设置指南	70

InPro 5000 (i) CO₂ Sensor Series
30 313 393 C

InPro 5000 (i) Senzory CO₂

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Obsah

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Quick Setup Guide

Contents

1	Introduction	15
2	Safety instructions	15
3	Product description	15
4	Sensor installation.....	16
4.1	Connection.....	16
5	Operation	16
5.1	Calibration of the pH sensor	16
5.2	Mounting the membrane body	17
5.3	Calibration of the CO ₂ sensor	17
6	Maintenance.....	17
7	Environmental protection.....	17
8	Warranty	17
9	Storage condition	17
10	EC Declaration of Conformity	17

InPro and ISM are trademarks of the METTLER TOLEDO Group.

1 Introduction

Thank you for buying this InPro™ 5000 (i) sensor from METTLER TOLEDO. InPro 5000 (i) sensors are intended solely for in-line measurement of dissolved CO₂ (carbon dioxide).

2 Safety instructions

 **Note:** The InPro 5000 (i) sensor should be installed, operated, maintained only by personnel familiar with the sensor and who are qualified for such work.

 **Note:** A defective sensor must neither be installed nor put into service.

 **Caution:** Only the maintenance work described in these operating instructions may be performed on the sensors.

 **Note:** When changing faulty components, use only original spare parts obtainable from your METTLER TOLEDO supplier.

 **Note:** No modifications to the sensor or the accessories are allowed. The manufacturer/supplier accepts no responsibility for damage caused by unauthorized modifications. The risk is borne entirely by the user.


3 Product description

Each sensor is supplied together with:

- 1 bottle CO₂-electrolyte 25 ml
- 1 quality control certificate
- 1 inspection certificate 3.1 (in accordance with EN 10204.3/3.1)
- 1 quick setup guide
- 1 EC Declaration of Conformity

Membrane body and cap sleeve are not mounted. The interior body is protected against drying out with a wetting cap, filled with CO₂-electrolyte.

Please check completeness.

1	VP connector	7	Retainer nut
2	ISM connector	8	Membrane body
3	Hexagonal screw	9	Interior body
4	Threaded sleeve Pg 13.5	10	Cap sleeve
5	Washer PTFE	11	Wetting cap
6	O-Ring (10.77 × 2.62 mm, Silicone FDA/USP VI)		

4 Sensor installation

4.1 Connection

The interior body (pH electrode) is ready for calibration.

- On unpacking, check the pH electrode for mechanical damage. Report any damage immediately to your METTLER TOLEDO supplier.
- Remove the watering cap, cap sleeve, membrane body and briefly rinse the sensor with de-ionized water.
- After rinsing, the sensor should only be dabbed dry with a tissue. Do not rub the pH-sensitive glass as this can lead to electrostatic charging and sluggish response times.
- Check inside the pH-sensitive glass for the presence of any air bubbles and remove same by gently swinging the electrode in a vertical plane (shaking down as with a thermometer).
- Calibrate the pH sensor as described in Section 5.
- Place the electrode into the housing as described in the instructions for the housing.

Note: Mounting upside-down or horizontally is not possible due to the internal pH electrode. An angle of 15° or greater to horizontal is needed.

- Connect the electrode and the pH transmitter with the corresponding cable (AK9 for the InPro 5000 (i), VP for the InPro 5000). Observe the connection scheme that accompanies the cable or the transmitter.
- ISM™ sensors such as InPro 5000i enable Plug and Measure startup and expanded diagnostics. For installation, commissioning and operation of the ISM system please refer to instructions for the transmitter, module and cable.

5 Operation

5.1 Calibration of the pH sensor

(Not needed for InPro 5000i)

a: Min. immersion depth > 25 mm / 1"

For the analog CO₂ sensor InPro 5000, we recommend a 2-point calibration in pH buffers 7.00 and 9.21. Before calibration, first remove wetting cap. Start with buffer pH 7.00 for determining zero point and then buffer pH 9.21 for determining the slope. Please refer to instruction manual of the transmitter for further details.

Note: On the transmitter, select Pt 1000 as the temperature sensor for the analog InPro 5000. The InPro 5000i does not require a temperature setting.

Note: Please do not leave the sensor in the calibration buffer for longer than 1 hour. Too long an exposure may result in high sensor drift in the process.

5.2 Mounting the membrane body

- After calibration of the pH electrode rinse the pH electrode with de-ionized water and carefully dab it with a paper tissue.
- Fill the new membrane body half way with CO₂ electrolyte.
- Make sure that all air bubbles are removed from the membrane body. Air bubbles can be removed by carefully tapping on the membrane body.
- Slip the membrane body over the pH electrode while holding the sensor in a vertical position. The excess electrolyte that will be displaced should be removed with a paper tissue.

Note: No electrolyte, sample media or contamination should be present between the membrane body and the cap sleeve. Please check carefully.

- Carefully slip the cap sleeve over the fitted membrane body, holding the sensor in a vertical position, and screw it tight. The cap sleeve must be clean and tight.

5.3 Calibration of the CO₂ sensor

The InPro 5000 is in situ sterilizable or in an autoclave. After sterilization a process calibration step is needed to achieve accurate measurements.

The process calibration can be performed by using the transmitter's calibration menu. For further details please refer to the main manual.

6 Maintenance

We recommend using a new membrane body after each batch. After each calibration of the pH electrode new CO₂ electrolyte should be filled in the membrane body.

If the pH electrode shows a sluggish response time and/or insufficient sensitivity, the pH-sensitive glass can be reactivated using an HF solution (p/n 51 319 053). For further details please refer to the main manual.

Typically, after 10–20 sterilization cycles the pH electrode should be replaced. Please refer to the instruction manual of the replacement interior body InPro 5000 (i).

7 Environmental protection

Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your Local Authority or retailer for recycling advice.

8 Warranty

Warranted free of manufacturing defects for a term of 12 months after delivery.

9 Storage condition

For storage we recommend filling the membrane body with fresh CO₂ electrolyte and reconnecting it to the sensor. The wetting cap, filled with CO₂-electrolyte, will prevent drying-out of the membrane.

10 EC Declaration of Conformity

The complete declaration is available in the declaration of conformity certificate.

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 Table composed in accordance with SJ/T 11364 (CN). This product is bearing the following symbol:

○：表示该有害物质在该部件所有均质材料中的含量均在GB/T 26572规定的限量要求以下
 ○：Indicates that the content of the hazardous substance in all homogeneous materials of the part is below the limit specified in GB / T 26572

x：表示该有害物质至少在该部件的某一均质材料中的含量超出GB/T 26572规定的限量要求。
 x：Indicates that the content of the toxic substance in at least one of the homogeneous materials of the part exceeds the limits specified in GB/T 26572.

* 仅限数字pH传感器 (InPro3XXXi)
 * For digital pH sensor only (InPro3XXXi)

For addresses of METTLER TOLEDO
Market Organizations please go to:
www.mt.com/contacts

METTLER TOLEDO Group
Process Analytics
Local contact: www.mt.com/contacts

Subject to technical changes
© 03/2023 METTLER TOLEDO
All rights reserved
Printed in Switzerland. 30 313 393 C

Management System
certified according to
ISO 9001 / ISO 14001

www.mt.com/pro

* 3 0 3 1 3 3 9 3 C *