



# CALIBRATION CERTIFICATE

Number VG 230210

Page 1 of 2

**Mettler-Toledo S.A.E.**  
Laboratorio de Calibración



## METTLER TOLEDO

Carrer Segrià 7 - 9  
08940 Cornellà de Llobregat (Barcelona)  
Tel. 93 223 76 00  
[Laboratorio.Calibracion@mt.com](mailto:Laboratorio.Calibracion@mt.com)

ITEM Piston operated burette

MANUFACTURER Mettler Toledo

MODEL DV 1005

NOMINAL VOLUME 5 ml

IDENTIFICATION B504536334

APPLICANT **METTLER-TOLEDO PAC RIM AG-TAIWAN**  
2F., No 17, Lane 171, Sec. 2 , Jiuzong Rd  
11494 Taipei City Taiwan

Calibration date 23 May 2023

Authorized Signatory

*This certificate in digital format is the original one.  
Any printing will be considered as a copy.*

Este certificado se expide de acuerdo con las condiciones de la acreditación concedida por ENAC, que ha comprobado las capacidades de medida del laboratorio y su trazabilidad a patrones nacionales o internacionales. ENAC es firmante del Acuerdo de Reconocimiento Mutuo (MLA) de calibración de European Cooperation for Accreditation (EA) y de International Laboratory Accreditation Cooperation (ILAC).  
*This certificate is issued in accordance with the conditions of accreditation granted by ENAC which has assessed the measurement capability of the laboratory and its traceability to national or international standards.*  
*ENAC is one of the signatories of the Multilateral Agreement of the European Cooperation for Accreditation (EA) and the International Laboratories Accreditation Cooperation (ILAC).*

**Instrument information**

Burette DV 1005  
 Serial number B504536334  
 Nominal volume 5 ml

**Calibration procedure**

Procedure PEC/MTE/22 based to the ISO 8655 norm and the METTLER TOLEDO manuals.  
 The measured volume corresponds with delivered volume (Ex) at the reference temperature of 20 °C.

**Calibration conditions**

Ambient temperature Min. 20,1 °C Max. 20,1 °C  
 Relative humidity 60,2 % Hr  
 Pressure 1010,9 mbar

Maintenance: Before calibration, seals (references 101003 and 25737) and piston have been replaced, and burette glass has been cleaned.

**Traceability**

Standard equipment used  
 Balance BAL01 AT201 - 5 decimal places balance  
 Burette drive VAL-T50 T50  
 Water temperature TER128 0,1 °C resolution  
 Ambient conditions REG02 (air Temp, rH)  
 Class III water 2307309810


The measured volume has been calculated conform to ISO/TR 20461:2023.  
 The traceability of measurements are referred to laboratories accredited by ILAC signatories accreditation bodies or national laboratories EURAMET signatories.

**Uncertainty**

The reported expanded uncertainty of measurement is stated as the Standard uncertainty of measurement multiplied by the coverage factor k = 2, which for a normal distribution corresponds to a coverage probability of approximately 95%. The standard uncertainty of measurement has been determined in accordance with document EA-4/02 M: 2022.

**Calibration results (after maintenance)**

Burette stroke	10%	30%	50%	100%	
<b>Selected volume</b>	<b>0,5</b>	<b>1,5</b>	<b>2,5</b>	<b>5</b>	<b>ml</b>
Water temperature	20,5	20,5	20,5	20,5	°C
Measured volume	499,2	1502,0	2502,9	5002,4	µl
<b>Systematic error</b>	<b>-0,8</b>	<b>2,0</b>	<b>2,9</b>	<b>2,4</b>	<b>µl</b>
Measurement uncertainty	2,5	2,5	2,5	2,5	µl
Random error	0,15	0,07	0,10	0,21	µl
Max. Permissible Sys. Error	15	15	15	15	µl
Max. Permissible Random Error	5	5	5	5	µl
<b>Evaluation*</b>	<b>Pass</b>	<b>Pass</b>	<b>Pass</b>	<b>Pass</b>	



\* Evaluation and permissible errors from ISO 8655-3:2022. According with this norm, uncertainty is not taken into account.

Mettler- Toledo S.A.E.  
 Carrer Segrià, 7-9  
 08940 Cornellà de Llobregat (Barcelona)

Technician: D. Gallardo  
 Remarks: ---